


PP3006 PRIMARY UTI CHROMOGENIC AGAR OPAQUE



Urinary Tract Infections (UTIs) account for 35-40% of all hospital acquired infections in the UK. Gram –ve aerobic bacteria are responsible for a significant proportion of UTIs with *E.coli* isolation rates at 80-90% of first time infections. Whilst this rate drops to approximately 70% in subsequent incidences of infections, a wide range of other pathogenic and opportunistic bacterial species can cause UTIs. Therefore, accurate diagnosis of the causative organism is required for the appropriate antibiotic to be prescribed.

Primary UTI Agar is a non-selective, chromogenic medium designed to facilitate the rapid detection and presumptive identification of the major pathogens responsible for UTIs. Based on the traditional CLED medium, to prevent the swarming of *Proteus* spp., two chromogenic compounds are included in the medium. The first chromogen, allows for the detection of *E.coli* (purple colonies) and the second chromogen allows for the detection of other coliforms in the KESC group of organisms (blue colonies). Tryptophan is also included in the medium to indicate the presence of Tryptophan Deaminase (TDA) activity allowing for the detection of *Proteus* spp, *Morganella* spp and *Providencia* spp.

This version of Primary UTI contains an additional agent in the chromogenic mix to provide an opaque, white background for improved colour contrast.

E.coli – Purple colonies

Klebsiella spp, *Enterobacter* spp, *Serratia* spp, *Citrobacter* spp. (KESC) – Metallic blue colonies

Staphylococcus spp. – White/pink colonies

Proteus spp, *Providencia* spp, *Morganella* spp. – Colourless colonies w brown halo

Enterococcus spp. – Turquoise blue colonies

Formula	gm/L	Properties	
Proprietary information	51.62	Appearance	Firm Gel
		Colour	White
		pH	6.7 ± 0.2
		Storage	2 - 8°C
		Shelf Life	70 days

Quality Control Test Organisms	Ref. No.	Result
<i>Enterococcus faecalis</i>	NCTC 12697	Turquoise / Blue Colonies
<i>Escherichia coli</i>	NCTC 12241	Purple Colonies
<i>Proteus mirabilis</i>	NCTC 10975	Clear Colonies with Brown Halo
<i>Enterobacter aerogenes</i>	NCTC 10006	Metallic Blue Colonies
<i>Staphylococcus aureus</i>	NCTC 12981	White Colonies

Recommended Incubation : Aerobically at 37°C ± 1°C for 18 - 24 hours