


EOrder.NET coming soon

Microsoft SQL Server 2005



Customer stock data is validated and stored in Database

ON-LINE
web access 24-7

ASP.NET / SOAP



E&O hosted server receives and validates data



Client computer

Stock Data is transferred from Barcode Scanner to client computer

Microsoft WCF



E&O's server stores customer data locally and sends to hosted server automatically




Client computer

Key Benefits

- Refined business logic with improved accuracy
- Upsized from Access to SQL 2008
- Automatic order importation on all E&O products
- Improved auditing including user/unit capture during transaction
- Remote hosting - accessed under licence from anywhere
- Operational tasks controlled by E&O - no backup requirements
- Integrated remote support - instant maintenance
- Lower support overhead for customer
- Lower client hardware requirement
- Simplified installation and up-upgrades, only one interface
- Customised reporting and charting options
- Flexible grouping, trusts, sites, departments
- Secure and password protected

Integrated Stock Management System



Overview

The third generation of E & O's integrated stock management system takes a significant technological leap onto the web. With a completely new web based front end, a secure, resilient hosted SQL 2008 database back end, and a much smaller, more robust client application to handle data logger transaction transfers.

Technologies

The system used the latest Microsoft .NET technologies and standards for each of the separate system components as listed below

- Microsoft SQL Server 2008
- Microsoft SQL Server Reporting Services
- ASP.NET / SOAP
- Microsoft WCF

